

# Polony Technology for Sequencing and Genotyping


Rob Mitra

Laboratory of George Church  
Harvard Medical School


## High Throughput Sequencing Applications

- De novo sequencing
- Resequencing (SNP discovery)
- SNP scoring
- Whole-genome mRNA abundance measurements
- Readout for genetic selections or screens


## Technology Overview: Polony Amplification



## Primer A has 5' Acrydite Modification




## Fluorescent In Situ Sequencing



## Advantages of Polony Sequencing

- Miniaturization (subpicoliter volumes)
- Parallel nature
- Integrated sample preparation
- Goal: 500 million bases per slide

Polony Radius can be as small as 5 microns





At this density, more than 5 million distinguishable polonies per slide



### Single Base Extensions on Gel-Immobilized Oligonucleotide Templates

Next Base	Added fl-dUTP	Fluorescent Intensity	Next Base	Added fl-dCTP	Fluorescent Intensity
<b>T</b>		100 +/- 6.0	<b>T</b>		0.55 +/- 0.1
<b>G</b>		0.59 +/- 0.5	<b>G</b>		0.42 +/- 0.1
<b>C</b>		0.58 +/- 0.4	<b>C</b>		100 +/- 5.3
<b>A</b>		0.18 +/- 0.2	<b>A</b>		0.30 +/- 0.1
<b>T</b>		1.5 +/- 1.0	<b>T</b>		3.4 +/- 0.5
<b>G</b>		100 +/- 6.2	<b>G</b>		5.6 +/- 0.8
<b>C</b>		0.8 +/- 0.3	<b>C</b>		8.8 +/- 0.8
<b>A</b>		2.4 +/- 0.4	<b>A</b>		100 +/- 10.8

## Incorporated Fluorescence Can be Removed Using a Modified Nucleotide




## Reaction Efficiency = 99.8%





## Sequencing on Colonies

7 8 9



ST30  
 3' P' TCACGAGT  
 5' P AGT**G**

ST30  
 3' P' TCACGAGT  
 5' P AGTG

ST30  
 3' P' TCACGAGT  
 5' P AGT**G****C**


Base Added: **G**

**T**

**C**

## Sequencing on Colonies

10 11 12



ST30  
 3' P' TCACGAGT  
 5' P AGT**G****C**

ST30  
 3' P' TCACGAGT  
 5' P AGT**G**

ST30  
 3' P' TCACGAGT  
 5' P AGT**G****C****T**

Base Added: **A**


**G**

**T**

# Sequencing on Colonies

13

14



ST30  
3'  P' TCACGAGT  
5'  P AGTGCTC

ST30  
3'  P' TCACGAGT  
5'  P AGTGCTCA

Base Added:

**C**

**A**

## Remaining Problems


- Gel attachment
- Calling multiple base runs (context dependence of fluorophore quantum yield)
- Cycle time
- Increase sequencing accuracy

## Other Applications of Polonies

- Genotyping single molecules
  - LOH analysis
  - Detecting rare SNPs (1 in  $10^3$ )
- Haplotyping
- Exon Analysis



## Direct Molecular Haplotyping with Polonies

- Determine phasing of SNPs over very long distances: Genotype (c/a t/g)




# Coamplify Two or More Loci and Genotype


Haplotype = ct  
ag





**Single molecule haplotypes (or RNA) detection:  
in situ amplified “polonies” (85% efficiency)**



## Exon Analysis



### Splice Variants



## Conclusions

- A large number of DNA template molecules can be amplified and manipulated in parallel using polony technology
- The amplified DNA can be extended by multiple rounds of nucleotide addition to give sequence information
- Polony technology may prove useful for genotyping, haplotyping, and exon typing

# Acknowledgements

George Church

Vincent Butty

Phil Busby

Aimee Dudley


The Church Lab

We are poster number 47,  
please stop by to see data from other projects




## Multiple polony DNA polymorphism linkage analysis

When two colored polonies overlap, the result is additive:

Red + Blue = Purple



# 1 Mb cis-polony DNA polymorphism linkage analysis






**Image Informatics**

1. Reject any similarly colored polonies closer together than 1/2 fully extended molecule length (eg, 165 μm ctr-to-ctr)
2. Reject any polony clusters that do not have exactly 4 members:
  - <4 = shear-damaged DNA or PCR artefact
  - >4 = random co-localization of template DNA molecules
3. Reject any clusters where polonies are not within predictable distances of each other

HMS Case 1929, Fig. 6

# Cis-polony mRNA “exon-typing” analysis



Reject any polony clusters that have only one member


**Results:**

- 01 B (yellow) reject
- 06 A-B (orange) 32%
- 05 B-C (green) 26%
- 05 A-C (purple) 26%
- 03 A-B-C (black) 16%

19 accepted 100%

HMS Case 1929, Fig. 17


## Polony Sequencing for SNP screening



15 base additions, 22 incorporation events


CYCLE: 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15  
BASE: A T A C A T G A C A C A G A C  
T

FITC (C) Channel



Mean Polony Intensity: 39.5

CY3 (T) Channel



Mean Polony Intensity: 6.5

Polony Template


3'  T A T T G T T A A A G T G T G T C C T T T G T C G A T A C T G G T A ...5'  
5'  A T A A C A A T T T C A C A C A G G A A A C

## 19 base additions, 26 incorporation events

CYCLE: 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19


BASE: A T A C A T G A C A C A G A C A G C <sup>C</sup><sub>T</sub>

FITC ( C ) Channel



Mean Polony Intensity: 0.25

CY3 ( T ) Channel




Mean Polony Intensity: 48.2


Polony Template

3' P T A T T G T T A A A G T G T C C T T T G T C G A T A C T G G T A ...5'  
5' P A T A A C A A T T T C A C A C A G G A A A C A G C T


## Polony Amplification



0



60



180



360

## Efficiency of Extension Reaction

Error accumulates exponentially with  
the number of incorporated bases

Assume reaction goes to 85% completion


After 6 extensions,  
 $(85\%)^6 = 38\%$  of the template molecules are extended correctly

The rest of the molecules are “out of phase”

## Future Direction

- Demonstrate sequencing out to 20 bases and determine error rate
- Automate sequencing protocol
- Demonstrate sequencing on colonies < 20 microns
- Demonstrate ability to haplotype alleles > 20kb apart

# Primer Extension 26 cycles, 34 Nucleotides



Mean Intensity: 58, 0.5

40, 6.5

0.3, 48


0.4, 43

Polony Template

3'  TATTGTTAAAGTGTGTCCTTTGTCGATACTGGTA...5'  
 5'  ATAACAATTTCACACAGGAAACAGCTATGACCAT



FITC ( C )

CY3 ( T )




# Proof of concept experiment

- Locus A polony (85% yield)
- Locus B polony (81% yield)
- Loci A+B cis-polony (69%)  
(polonies are ~300  $\mu$  diam)




**Double Cut**  
(EcoRI & NcoI)



**Single Cut**  
(EcoRI only)

HMS Case 1929, Fig. 3


## Single molecule haplotypes (or RNA) detection: in situ amplified “polonies” (85% efficiency)



Extension estimated to go to 99.8% completion

Extension with unlabeled dTTP	Next Base	Polymerase Trapping	Fluorescent Intensity
0 seconds	T		100 +/- 14.5
	C		1.83 +/- 0.024
60 seconds	T		1.97 +/- 1.38
	C		1.73 +/- 1.22

## Technology Overview: Polony Sequencing



# High Throughput Sequencing

Millions of DNA fragments



Clone



Amplify



Sequence



Processed